



Position: Quality Assurance Manager
Reports To: Technical Director
Location: Castellar del Vallès

Company description:

Since its founding in Chicago (USA) more than 100 years ago, ITW has become one of the world's leading diversified manufacturers of specialized industrial equipment, consumables and related service businesses.

ITW business serve local customers and markets around the globe. The company has operations in 56 countries that employ more than 50.000 people who adhere to the highest ethical standards. ITW has 85 businesses divided into 7 Segments, which includes: Welding, Automotive OEM, Food Equipment Group, Construction Products, Polymers & Fluids, Test and Measurement & Electronics, and Specialty.

Summary

As the Quality Assurance Manager for Castellar del Vallès you will be responsible to ensure that all external and internal requirements are met before our product reaches our customers. You will be responsible for inspecting procedures and outputs and identifying mistakes or non-conformity issues. You're responsible of the definition, implementation and follow-up, of the integrated management system.

What you will do

- Maintaining all documents about assuring the quality in products and services established and updated.
- Assure the quality of the products and services that the company gives to its clients.
- Comply with procedures and quality rules and the ones determined as general character.
- Coordinate and support on-site audits conducted by external providers.
- Interpret and implement quality assurance standards.
- Devise sampling procedures and direction for recording and reporting quality.
- Investigate customer complaints and non-conformance issues.
- Identify training needs and organize training interventions to meet quality standards.
- Responsible for document management systems.
- Assure ongoing compliance with quality and industry regulatory requirements.

Division Description:

The ITW Reagents Division has been established by the companies Panreac Química SLU (Spain), AppliChem GmbH (Germany) and Nova Chimica (Italy) with the aim of bringing together the experience and know-how of three seasoned companies.

The Division offers the most innovative and best quality products for chemical, biological, pharmaceutical and medical research and production applications.

Nowadays ITW Reagents markets its products worldwide through an extensive distribution network in more than 80 countries under the brand PanReac AppliChem.

What are we looking for

- Proven knowledge and understanding of the regulatory Pharma market and experience of the processes of Quality Assurance
- Excellent interpersonal, communication, and presentation skills
- Strong goal orientation, with 80/20 mindset
- Ability to persuade and convince others, according to the Company needs, objectives and policies.
- Attention to detail, communication skills, data collection.
- Problem analysis and problem solving.

Education and Experience

- Bachelor's degree in Chemistry, chemical engineer or similar, certification an advantage including Quality Auditor, Quality Engineer, Quality Improvement Associate, Six Sigma.
- Minimum 2 years experience as a Quality Assurance Manager in Pharma companies.
- Experience in a Quality Assurance Dpt. In a chemical and/or food company is a plus.
- Fluent in Spanish & English (spoken and written). German is a plus.
- Knowledge of tools, concepts and methodologies of QA. Knowledge of relevant regulatory requirements.

Contact details:

HR Department

E-mail: hr.es@itwreagents.com